LAB3_Tips

SURVIVAL TIPS

File formats (really important!!)

Bioinformatics use different file formats and each of these formats is characterized by specific features and fields. Let's see someone.

FASTA

A sequence record in a FASTA format consists of a single-line description (sequence name), followed by line(s) of sequence data. The first character of the description line is a greater-than (">") symbol.

>read id 0

 ${\tt GGTATGCTTCTGGGGCGGCAGTCGATAGGGCTAGACTCAGGTCCCGTGGC}$

>read id 1

CACTGTGGCCCTCTTGGGGGGTGTCCACACGCCGCCCGTCGGCCCCCTCC

>read id 2

GTTCTGTGGGTACCTCGCGGTTATGGTGTCGGGGGTATCCAAGGCACCCC

FASTQ

Similar to FASTA file, but with mapping quality information for each base. Both the sequence letter and quality score are each encoded with a single ASCII character.

@ERX288614.1 HWI-ST1362:33:D1J0JACXX:6:1101:1687:2354 length=101 TTTTTCTAGACGGCAGGTCAGGTCCACCACTGACACGTTGGCAGTGGGACACGGAAGGCCATGCCAGTGAGCTTCCCGTTC AGCTCAGGGATGACCTTGC

SAM

SAM stands for Sequence Alignment/Map format and it is the most common file output format for aligners. It is a TAB-delimited text format consisting of a header section, which is optional, and an alignment section.

@SQ			LN:1337	197422			
@SQ			SN:18		LN:80373285		
@PG	ID:bwa		PN:bwa	VN:0.7.17-r1	188 CL:bwa	mem	
/home/marta/Documents/BIOINFORMATICS/BioInfoCourse/LAB_alignment/tools/bwa_index/b							
wa_index			mate_1.fq		mate	_2.fq	
ERX288614.1	99	10	55667614	39	12S79M1I9M	=	

55667850 318

TTTTTCTAGACGGCAGGTCAGGTCCACCACTGACACGTTGGCAGTGGGGACACGGAAGGCCATGCCAGTGAGCTTCCCGTTC AGCTCAGGGATGACCTTGC

-318

AGTCCTTCCACGATACCAAAGTTGTCATGGATGTCCTTGGCCAGGGGTGCTAAGCAGTTGGTGCAGGAGGCATTGCTGA
TGATCTTGAGGCTGTTGTC

Remember that an aligner can report multiple alignments for the same read!! Depending from the application, this could be an issue to be fixed.

Take a look here https://samtools.github.io/hts-specs/SAMv1.pdf for all details about this file format, paying particular attention to section 1.4 about mandatory fields (page 6).

Col	Field	Type	Regexp/Range	Brief description
1	QNAME	String	[!-?A-~]{1,254}	Query template NAME
2	FLAG	Int	$[0, 2^{16} - 1]$	bitwise FLAG
3	RNAME	String	* [:rname:^*=][:rname:]*	Reference sequence NAME ⁹
4	POS	Int	$[0, 2^{31} - 1]$	1-based leftmost mapping POSition
5	MAPQ	Int	$[0, 2^8 - 1]$	MAPping Quality
6	CIGAR	String	* ([0-9]+[MIDNSHPX=])+	CIGAR string
7	RNEXT	String	* = [:rname:^*=][:rname:]*	Reference name of the mate/next read
8	PNEXT	Int	$[0, 2^{31} - 1]$	Position of the mate/next read
9	TLEN	Int	$[-2^{31}+1, 2^{31}-1]$	observed Template LENgth
10	SEQ	String	* [A-Za-z=.]+	segment SEQuence
11	QUAL	String	[!-~]+	ASCII of Phred-scaled base QUALity+33

N.B. Genome positions in SAM files are in **1-based coordinate system**

BAM

BAM is the binary version of a SAM file. This means that BAM files are smaller than SAM files and this property is really helpful especially when we are working with huge files (e.g. a 30 GB SAM file can be compressed into a 17 GB BAM file). However, since BAM files are binary files, they are not human readable.

N.B. Genome positions in BAM files are in **0-based coordinate system.**

GTF

The Gene transfer format (GTF) is a file format used to hold information about gene structure of a reference genome. It is a tab-delimited text format based on the general feature format (GFF). This file format is really helpful when you want to know which biological feature (gene, exon,CDS, ...) is present in which genome positions.

```
#!genebuild-last-updated
                                                                               2016-06
                                      14409
          havana
                   gene
                             11869
                                                                               gene_id
"ENSG00000223972"; gene_version "5"; gene_name "DDX11L1"; gene_source "havana";
gene_biotype
                        "transcribed_unprocessed_pseudogene";
                                                                           havana gene
"OTTHUMG00000000961";
                                          havana gene version
                                                                                   "2":
         havana transcript
                                   11869
                                            14409
                                                                               gene_id
1
                                         "5";
                                                                    "ENST00000456328";
"ENSG00000223972";
                       gene version
                                                 transcript id
transcript_version "2"; gene_name "DDX11L1"; gene_source "havana"; gene_biotype
"transcribed_unprocessed_pseudogene";
                                                                 "OTTHUMG00000000961";
                                              havana_gene
havana_gene_version "2"; transcript_name "DDX11L1-002"; transcript_source "havana";
transcript biotype "processed transcript"; havana transcript "OTTHUMT00000362751";
                                           "basic";
havana_transcript_version "1";
                                    tag
                                                     transcript_support_level
    havana exon 11869 12227 .
                                    +
                                            gene_id "ENSG00000223972"; gene_version "5";
1
transcript_id "ENST00000456328"; transcript_version "2"; exon_number "1"; gene_name "DDX11L1";
gene source "havana"; gene biotype "transcribed unprocessed pseudogene"; havana gene
"OTTHUMG00000000961";
                                                 "2";
                                                                         "DDX11L1-002";
                          havana gene version
                                                       transcript name
transcript source
                  "havana";
                             transcript_biotype
                                                "processed transcript";
                                                                       havana transcript
"OTTHUMT00000362751";
                         havana_transcript_version
                                                    "1";
                                                          exon id
                                                                   "ENSE00002234944":
                   "1":
                                        "basic":
exon_version
                              tag
                                                      transcript_support_level
    havana exon
                 12613 12721 .
                                   +
                                        .
                                            gene_id "ENSG00000223972"; gene_version "5";
transcript_id "ENST00000456328"; transcript_version "2"; exon_number "2"; gene_name "DDX11L1";
gene_source "havana"; gene_biotype "transcribed_unprocessed_pseudogene"; havana_gene
"OTTHUMG00000000961";
                                                "2";
                          havana gene version
                                                       transcript name
                                                                         "DDX11L1-002";
transcript source
                  "havana":
                             transcript biotype
                                                "processed transcript";
                                                                       havana transcript
"OTTHUMT00000362751";
                          havana_transcript_version
                                                    "1";
                                                          exon id
                                                                    "ENSE00003582793";
exon_version
                   "1";
                              tag
                                        "basic";
                                                      transcript_support_level
                                                                                    "1";
```

N.B. Genome positions in GTF files are in 1-based coordinate system

VCF

VCF is a text file format (most likely stored in a compressed manner) and **it holds information about reads**. It contains meta-information lines, a header line, and then data lines each containing information about a position in the genome. The format also has the ability to contain genotype information on samples for each position.

```
##fileformat=VCFv4.2
##FILTER=<ID=PASS,Description="All filters passed">
##bcftoolsVersion=1.9+htslib-1.9
##bcftoolsCommand=mpileup --fasta-ref reference chr10 chr18.fa sorted.bam
##reference=file://reference chr10 chr18.fa
##contig=<ID=10,length=133797422>
##contig=<ID=18,length=80373285>
##ALT=<ID=*, Description="Represents allele(s) other than observed.">
##INFO=<ID=I16, Number=16, Type=Float, Description="Auxiliary tag used for calling,
see description of bcf callret1 t in bam2bcf.h">
##INFO=<ID=QS, Number=R, Type=Float, Description="Auxiliary tag used for calling">
##FORMAT=<ID=PL,Number=G,Type=Integer,Description="List of Phred-scaled genotype
likelihoods">
#CHROM
                                                        INFO
          POS
                 ID
                                     QUAL
                                              FILTER
                                                                FORMAT
                                                                          sorted.bam
<mark>10</mark>
              48636
                                           Т
                                                                       0
DP=1;I16=0,0,0,0,0,0,0,0,0,0,0,0,0,0,0,0;QS=0,0;MQ0F=0
                                                                  0,0,0
```

```
10
             48637
                                                                     0
DP=2;I16=0,0,0,0,0,0,0,0,0,0,0,0,0,0,0;QS=0,0;MQ0F=0
                                                                0,0,0
             48638
                                                                     0
DP=2;I16=1,0,0,0,70,4900,0,60,3600,0,0,2,4,0,0;QS=1,0;MQ0F=0
                                                                  PL
                                                                        0,3,60
10
DP=2;I16=1,0,0,0,74,5476,0,0,60,3600,0,0,3,9,0,0;QS=1,0;MQ0F=0
                                                                  PL
                                                                         0,3,60
                                                                     0
DP=2;I16=1,0,0,0,70,4900,0,0,60,3600,0,0,4,16,0,0;QS=1,0;MQ0F=0
                                                                   PL
                                                                         0,3,60
             48641
                                                                     0
DP=2;I16=1,0,0,0,70,4900,0,0,60,3600,0,0,5,25,0,0;QS=1,0;MQ0F=0
                                                                   PL
                                                                         0,3,60
                                                     T,<*>
DP=2;I16=0,0,1,0,0,0,70,4900,0,0,60,3600,0,0,6,36;QS=0,1,0;SGB=-0.379885;MQ0F=0
PL
      60,3,0,60,3,60
```

INSTALL SAMTOOLS, BCFTOOLS and download human gtf file

```
conda activate Bioinfo_labs (or source activate Bioinfo_labs)
conda install -c bioconda samtools
conda install -c bioconda bcftools
```

SAMTOOLS, really basic usage

SAM Tools provide various utilities for manipulating alignments in the SAM format, including sorting, merging, indexing and generating alignments in a per-position format.

SAM/BAM conversions

As previously pointed out, BAM format is the binary version of a SAM file. The conversion can be performed using samtools:

```
samtools view -S -b my.sam > my.bam
```

BAM sorting

When multiple selections have to be performed onto a huge file, it is convenient to sort that file according to certain criteria (e.g. genomic region) in order to search for the required information in a faster way. For BAM sorting you can use:

```
samtools sort my.bam > my-sorted.bam
```

FLAGS in a SAM/BAM file

Flags are used to keep track of alignment information in a compact way and uses 12 bits. Let's see how they are computed.

FLAG: Combination of bitwise FLAGs. Each bit is explained in the following table:

Bit position			
from rigth	Bit		Description
0	1 0x1		template having multiple segments in sequencing
1	2	0x2	each segment properly aligned according to the aligner
2	4	0x4	segment unmapped
3	8	0x8	next segment in the template unmapped
4	16	0x10	SEQ being reverse complemented
5	32	0x20	SEQ of the next segment in the template being reverse complemented
6	64	0x40	the first segment in the template
7	128	0x80	the last segment in the template
8	256	0x100	secondary alignment
9	512	0x200	not passing filters, such as platform/vendor quality controls
10	1024	0x400	PCR or optical duplicate
11	2048	0x800	supplementary alignment

Some examples with FLAGS:

00000000001 --> $2^0 = 1$ --> template having multiple segments in sequencing

00000000010 \rightarrow 2¹ = 2 \rightarrow each segment properly aligned according to the aligner

 $00000000100 --> 2^2 = 4 -->$ segment unmapped

 $0001000000000 --> 2^8 = 256 --> secondary alignment$

Bits can be combined:

00000001100 --> $2^2 + 2^3 = 4 + 8 = 12$ --> segment unmapped **and** next segment in the template unmapped

Using samtools to filter alignments in a SAM/BAM file using FLAGS

Among the many potentials of samtools view there is that of filtering the reads using the FLAG field in SAM files.

To have an overall idea of how samtools view works open your terminal, activate Bioinfo_labs environment and type

samtools view

to get its manual page. Take a look at -f and -F options. These options allow us to filter reads in a SAM file following a combination of criteria relying on alignment flags.

- -f INT only include reads with all of the FLAGs in INT present
- -F INT only include reads with none of the FLAGS in INT present

E.g. to obtain a SAM file with no unmapped reads and no secondary alignments we can exploit bit in position 2 (2^2 =4) and bit in position 8 (2^8 =256).

samtools view -F 260 bwa_out.sam > unique_aligned.sam

-F 260 $(2^{\frac{2}{3}} + 2^{\frac{8}{3}})$ means that in the final SAM file will be printed only reads for which bit number 2 or bit number 8 is not set to 1.

BCFTOOLS

BCFtools is a set of utilities that manipulate variant calls in the Variant Call Format (VCF) and its binary counterpart BCF.

BCFTOOLS to create VCF file

Use bcftool to convert a sorted BAM file into VCF:

bcftools mpileup --fasta-ref reference_chr10_chr18.fa sorted.bam >
sorted.vcf